

## Tacta<sup>®</sup> Pipette Multipacks The Smart Solution for Your Lab

### Tacta<sup>®</sup> Mechanical Pipettes

- Reliable results, even over lengthy pipetting periods
- Comfortable to hold due to the ergonomically designed handle
- Low pipetting and tip ejection forces
- Controlled and smooth tip ejection with the new Sartorius Optiject<sup>®</sup> technology
- The unique Sartorius Optilock<sup>®</sup> system provides flexibility for volume adjustment and locking
- Calibration adjustment to provide accurate results for various liquid types
- Easy to clean; with only three parts to disassemble
- Fully autoclavable

### Sartorius Optifit Tips

- High-quality standard tips
- Excellent solution for various applications
- Versatile selection of packages and purity options



# Tacta® Multipacks

Multipack Order Code	Tacta® Pipettes	Optifit Tips	Accessories
Tacta® Pipette 3-pack 10 LH-729670	1-channel <ul style="list-style-type: none"> <li>▪ 0,5 – 10 µL</li> <li>▪ 10 – 100 µL</li> <li>▪ 100 – 1,000 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,1 – 10 µL</li> <li>▪ 0,5 – 200 µL</li> <li>▪ 10 – 1,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 3-pack 20 LH-729671	1-channel <ul style="list-style-type: none"> <li>▪ 2 – 20 µL</li> <li>▪ 20 – 200 µL</li> <li>▪ 100 – 1,000 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,5 – 200 µL</li> <li>▪ 10 – 1,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 4-pack 20 LH-729672	1-channel <ul style="list-style-type: none"> <li>▪ 0,5 – 10 µL</li> <li>▪ 2 – 20 µL</li> <li>▪ 20 – 200 µL</li> <li>▪ 100 – 1,000 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,1 – 10 µL</li> <li>▪ 0,5 – 200 µL</li> <li>▪ 10 – 1,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 4-pack 100 LH-729673	1-channel <ul style="list-style-type: none"> <li>▪ 0,5 – 10 µL</li> <li>▪ 10 – 100 µL</li> <li>▪ 20 – 200 µL</li> <li>▪ 100 – 1,000 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,1 – 10 µL</li> <li>▪ 0,5 – 200 µL</li> <li>▪ 10 – 1,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 5-pack 10 LH-729674	1-channel <ul style="list-style-type: none"> <li>▪ 0,5 – 10 µL</li> <li>▪ 10 – 100 µL</li> <li>▪ 20 – 200 µL</li> <li>▪ 100 – 1,000 µL</li> <li>▪ 500 – 5,000 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,1 – 10 µL</li> <li>▪ 0,5 – 200 µL</li> <li>▪ 10 – 1,000 µL</li> <li>▪ Tip Rack – 50 tips</li> <li>▪ 100 – 5,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 5-pack 20 LH-729675	1-channel <ul style="list-style-type: none"> <li>▪ 2 – 20 µL</li> <li>▪ 10 – 100 µL</li> <li>▪ 20 – 200 µL</li> <li>▪ 100 – 1,000 µL</li> <li>▪ 500 – 5,000 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,5 – 200 µL</li> <li>▪ 10 – 1,000 µL</li> </ul> Tip Rack – 50 tips <ul style="list-style-type: none"> <li>▪ 100 – 5,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 3+1-pack LH-729676	1-channel <ul style="list-style-type: none"> <li>▪ 0,5 – 10 µL</li> <li>▪ 20 – 200 µL</li> <li>▪ 100 – 1,000 µL</li> </ul> 8-channel <ul style="list-style-type: none"> <li>▪ 30 – 300 µL</li> </ul>	Tip Racks – 96 tips <ul style="list-style-type: none"> <li>▪ 0,1 – 10 µL</li> <li>▪ 0,5 – 200 µL</li> <li>▪ 5 – 350 µL</li> <li>▪ 10 – 1,000 µL</li> </ul>	Linear Stand
Tacta® Pipette 3-pack max LH-729677	1-channel <ul style="list-style-type: none"> <li>▪ 100 – 1,000 µL</li> <li>▪ 500 – 5,000 µL</li> <li>▪ 1,000 – 10,000 µL</li> </ul>	Tip Racks <ul style="list-style-type: none"> <li>▪ 10 – 1,000 µL – 96 tips</li> <li>▪ 100 – 5,000 µL – 50 tips</li> <li>▪ 100 – 10,000 µL – 35 tips</li> </ul>	Linear Stand



## Germany

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen  
Phone +49 551 308 0

For further information, visit  
[www.sartorius.com](http://www.sartorius.com)

## USA

Sartorius Corporation  
565 Johnson Avenue  
Bohemia, NY 11716  
Phone +1 631 254 4249  
Toll-free +1 800 635 2906

## Finland & Baltics

Sartorius  
Biohit Liquid Handling Oy  
Laippatie 1  
00880 Helsinki  
Phone +358.9.755.951  
Fax +358.9.755.95.292  
[linfo.finland@sartorius.com](mailto:linfo.finland@sartorius.com)